

Data and text mining

***Lipid-Pro*: a computational lipid identification solution for untargeted lipidomics on data-independent acquisition tandem mass spectrometry platforms**

Zeeshan Ahmed^{1,2,*}, Michel Mayr³, Saman Zeeshan², Thomas Dandekar², Martin J. Mueller³ and Agnes Fekete³

¹Department of Neurobiology and Genetics, Biocenter, ²Department of Bioinformatics, Biocenter and ³Department of Pharmaceutical Biology, Biocenter, University of Wuerzburg, Wuerzburg, Germany

*To whom correspondence should be addressed.

Associate Editor: Jonathan Wren

Received on July 15, 2014; revised on November 24, 2014; accepted on November 25, 2014

Abstract

Summary: A major challenge for mass spectrometric-based lipidomics, aiming at describing all lipid species in a biological sample, lies in the computational and bioinformatic processing of the large amount of data that arises after data acquisition. *Lipid-Pro* is a software tool that supports the identification of lipids by interpreting large datasets generated by liquid chromatography–tandem mass spectrometry (LC–MS/MS) using the advanced data-independent acquisition mode MS^E. In the MS^E mode, the instrument fragments all molecular ions generated from a sample and records time-resolved molecular ion data as well as fragment ion data for every detectable molecular ion. *Lipid-Pro* matches the retention time-aligned mass-to-charge ratio data of molecular- and fragment ions with a lipid database and generates a report on all identified lipid species. For generation of the lipid database, *Lipid-Pro* provides a module for construction of lipid species and their fragments using a flexible building block approach. Hence, *Lipid-Pro* is an easy to use analysis tool to interpret complex MS^E lipidomics data and also offers a module to generate a user-specific lipid database.

Availability and implementation: *Lipid-Pro* is freely available at: <http://www.neurogenetics.biocenter.uni-wuerzburg.de/en/project/services/lipidpro/>

Contact: zeeshan.ahmed@uni-wuerzburg.de

Supplementary information: [Supplementary data](#) are available at *Bioinformatics* online.

1 Introduction

Lipids are compounds of biological origin that are soluble in non-polar solvents. In biological systems, lipids perform ubiquitous functions such as construction of cell membranes, stabilization of membrane-bound proteins, balancing of energy metabolism or cell signaling (Wymann and Schneider, 2008). All lipid species in a cell or tissue are referred to as the lipidome. Changes within lipidomes have been closely linked to stress responses or various diseases

including diabetes, obesity, heart diseases or neurodegenerative diseases (reviewed in Murphy and Nicolaou, 2013).

Liquid-chromatography coupled to tandem mass spectrometry (LC–MS/MS) is suitable for lipidome analysis since it can cover a broad range of lipids (Harald *et al.*, 2012). Both, data-dependent acquisition (DDA) and data-independent acquisition (DIA) are widely used for LC–MS-based lipidomics. In DDA, the most abundant molecular ions at each retention time (rt) are successively




Fig. 1. The graphical user interfaces of the five principal modules (*Building Block*, *Precursor Candidates*, *Fragment Candidates*, *Mass Analyzer* and *Lipid Identifier*) of *Lipid-Pro* are shown. Results can be exported in Microsoft Excel format

selected for isolation and fragmentation. Despite the advantages of the DDA approach, there are some inherent limitations such as often irreproducible molecular ion selection, undersampling and long instrument cycle times. DIA overcomes those limitations. MS^E, supported by Waters qTOF instrument platforms, acquires intact molecular- and fragment ions in an unbiased manner by continuously switching the energy for collision-induced dissociation between low (MS) and high (MS/MS) energy. DIA analysis requires coherent and intricate data processing since it produces complex data containing information on all molecular and fragment ions at each rt.

Commercial and free downloadable software such as mzMine 2 (Pluskal *et al.*, 2010) or Lipid Data Analyser (Hartler *et al.*, 2011) provide support to align and compare the spectral features of LC–MS data. These applications identify lipid species in LC–MS raw data by matching their monoisotopic masses or elemental compositions with a custom-based or online database. However, these software applications do not consider fragmentation data. LipidBlast (Kind *et al.*, 2013) identifies lipid species by matching the measured MS/MS spectra to an *in silico* generated library after converting raw data into the MGF format. LipidBlast can be used for the analysis of MS^E data, but it is laborious to match hundreds of different lipid species of the database with LC–MS/MS data.

We developed the *Lipid-Pro* software solution for rapid identification of lipids in preprocessed DIA data acquired with LC–MS/MS using rt-aligned monoisotopic masses of molecular ions (MS) and fragment ions (MS/MS). Since *Lipid-Pro* provides a module for construction of lipid species and their characteristic fragments using a flexible building block approach, rare or even non-described lipids beside the common species can be identified in total lipid extracts

of biological samples. To meet the technological objectives of this research, we took a step forward in the development of a new user friendly, modular and client-based database management system. The *Lipid-Pro* application can be configured on Microsoft Windows platforms following a simple six-step installation process.

2 Lipid-Pro

We developed *Lipid-Pro* as a desktop application (Fig. 1), helpful for the identification of pre-defined lipid species (termed candidates) in large LC–MS^E datasets (or other DIA datasets from other instrument platforms provided in the *.csv data format). For automatic generation of a lipid database, *Lipid-Pro* provides modules for the construction of complex lipid species and their fragments using a building block approach (Yang *et al.*, 2009). By comparison with the self-generated lipid database, pre-defined lipids are identified in the preprocessed datasets by matching the rt-aligned mass-to-charge ratios (*m/z*) of lipid molecular ions and fragment ions in an efficient manner. As output, *Lipid-Pro* delivers information on the identified lipid species including name, rt, *m/z* as well as the peak areas of lipid species and their fragments.

Lipid-Pro comprises five modules working in product-line architecture: *Building Block*, *Precursor Candidates*, *Fragment Candidates*, *Mass Analyzer* and *Lipid Identifier* (Fig. 2). Each module performs its task independently and its output is used as an input for the next module until lipid identification is achieved.

2.1 Lipid database

The majority of lipid molecular species consist of a linear combination of a small number of chemical entities or building blocks from

- Murphy, S.A. and Nicolaou, A. (2013) Lipidomics applications in health, disease and nutrition research. *Mol. Nutr. Food Res.*, **57**, 1336–1346.
- Pluskal, T. *et al.* (2010) MZmine2: modular framework for processing, visualizing, and analyzing mass spectrometry-based molecular profile data. *BMC Bioinformatics* **11**, 395–405.
- Yang, K. *et al.* (2009) Automated lipid identification and quantification by multidimensional mass spectrometry-based shotgun lipidomics. *Anal. Chem.*, **81**, 4356–4368.
- Wymann, M.P. and Schneider, R. (2008) Lipid signalling in disease. *Nat. Rev. Mol. Cell Biol.*, **9**, 162–176.